




# SPLINX RETAINERS



Permanent Retention made Easy!

- The first bonded lingual retainer wire that does not fray when cut.
- Made from a flattened coaxial wire to a precise dimension of .20 mm x .75 mm
- Lays neatly on the lingual surface of the teeth for precise bonding and retention.
- Works well with all types of adhesive and is easily formed to fit the tooth surface.
  - Suitable for Maxillary and Mandibular arches.



SPLINX™ PACK OF EIGHT 8 INCH STRIPS  
ITEM# 550-120 \$61/EACH



Contact Your Lancer Representative For Further Information

1-800-854-2896

Website: [www.lancerortho.com](http://www.lancerortho.com)